

BONDING & BANDING



Straight Bracket Remover- Anterior

Quality construction and design allows easy removal of bonded brackets without patient discomfort. Precision diamond honed edges work on all types of brackets. Tip width is 3 mm.



Angulated Bracket Remover- Wide

Conveniently angled head for working from behind the patient when removing brackets. Precision edges allow for easy removal of brackets and excess bonding material. 6 mm double-wide inserted jaws.



Posterior Band Remover

Our most popular size band remover has an inserted tip which designed to wedge under the band while the 5 mm replaceable Teflon bumper rests on the occlusal surface for easy band removal.



Crown & Band Crimping Plier

Slim tip design makes it simple to crimp and reform crowns and bands for a snug fit. Great control with this handy instrument.



Torquing Set

These pliers replace torquing keys and allow you to apply torque without distorting other segments of the archwires. Non-slip grip tips allow you to place labial or lingual torque on adjoining teeth. Available as a set or individually.

Wire Range:

.022" x .028" Rectangular

LAB PLIERS



Adams Plier

Two .040" square beaks with parallel opening at 1 mm. Made on our heavy duty lab forging with slightly larger tip size.



Hook Crimping Plier

The head of this plier is designed with a slot to fit over the base of the hook for a secure hold. One squeeze and the hook is secured on the archwire.



Crimp Plier/Angled

Angled head is designed to improve access to the posterior region. One squeeze and the hook is secure.
